

Box gross weight 18,4 kg – Vol. box 0,15 m³

Mesh back operative chair

Structure: 12mm thick plywood inner seat. 40 kg/m³ density fire retardant injected foam seat padding. Nylon backrest upholstered in RV mesh with height adjustable plastic lumbar support.

Armrests: **BR020.N** - black nylon height adjustable armrests with PU armpad, 90 mm stroke.



BR021.N - black nylon 3D adjustable armrests with PU armpad, 90 mm stroke.

Mechanism: **MS11A** – synchronized mechanism with 4 locking positions and side tension control. With seat slide 70 mm stroke.



MS1AE – synchronized mechanism with 4 locking positions and side tension control.

Base: **BG018.N** – Ø 700 mm black nylon base with metal ring.



BG004.L – Ø 700 mm polished aluminum base.



Castors: **R26** – Ø 65 mm self-breaking castors for soft floors.




R27 – Ø 65 mm self-breaking castors for hard floors.




Seat slide adjustment.

Pull the lever on the left hand side of the mechanism to adjust the seat depth.



Seat height adjustment.

Pull the right hand side lever while seated to lower the seat height. Pull the lever while standing to increase the seat height.



Tilt adjustment.

Twist the left hand side knob to release the synchronized action. Twist the knob back to lock the synchronized in the desired position.



Tilt tension adjustment.

Twist the right hand side knob to adjust the tilt tension.



Lumbar support adjustment.

Lift or lower the lumbar support in the desired position.

