

Box gross weight 18,4 kg – Vol. box 0,15 m³

Padded back operative chair

Structure:

12mm thick plywood inner seat. 40 kg/m³ density fire retardant injected foam seat padding. Nylon upholstered backrest padded in polyurethane foam with height adjustable plastic lumbar support.

Armrests:



BR020.N - black nylon height adjustable armrests with PU armpad, 90 mm stroke.

BR021.N - black nylon 3D adjustable armrests with PU armpad, 90 mm stroke.

Mechanism:



MS11A – synchronized mechanism with 4 locking positions and side tension control. With seat slide 70 mm stroke.

MS1AE – synchronized mechanism with 4 locking positions and side tension control.

Base:



BG018.N – Ø 700 mm black nylon base with metal ring.



BG004.L – Ø 700 mm polished aluminum base.

Castors:




R26 – Ø 65 mm self-breaking castors for soft floors.



R27 – Ø 65 mm self-breaking castors for hard floors.

Seat slide adjustment.

Pull the lever on the left had side of the mechanism to adjust the seat depth.

**Seat height adjustment.**

Pull the right hand side lever while seated to lower the seat height. Pull the lever while standing to increase the seat height.

**Tilt adjustment.**

Twist the left hand side knob to release the synchronized action. Twist the know back to lock the synchronized in the desired position.

**Tilt tension adjustment.**

Twist the right had side knob to adjust the tilt tension.

**Lumbar support adjustment.**

Lift or lower the lumbar support in the desired position.

