

Box gross weight 17,1 kg – Vol. box 0,15 m³

Operative chair

Structure: Plastic inner seat with metal bars. 80 mm stroke height adjustable backrest. 30 kg/m³ density fire retardant polyurethane foam padding.

Armrests: BR020.N - black nylon height adjustable armrests with PU armpad, 90 mm stroke.



BR021.N - black nylon 3D adjustable armrests with PU armpad, 90 mm stroke.



BR185.N - fixed polypropylene armrests.

Mechanism: MS11H – Synchronized mechanism with 4 locking positions and knob for tilt tension adjustment.



Base: BG055.N – Ø 650 mm 5 star nylon base with metal ring.



Castors: R04 – Ø 50 mm self-breaking castors for soft floors.




R05 – Ø 50 mm self-breaking castors for hard floors.



Seat height adjutment.

Pull the right hand side lever while seated to lower the seat height. Pull the lever while standing to increase the seat height.

**Tilt adjutment.**

Release the left hand side lever to unlock the tilt action. Lock the lever to lock the tilt action.

**Tilt tension adjutment.**

Twist the knob clockwise to increase the tilt tension. Twist it anti-clockwise to decrease the tilt tension.

**Backrest height adjustment.**

Pull the backrest upwards to increase the backrest height. Pull the backrest at its maximum height to set the backrest at its lowest position.

