

Box gross weight 19,2 kg – Vol. box 0,25 m³

Operative chair

Structure:

Plastic inner seat with metal bars, 60 kg/m³ density fire retardant injected foam padding 50 mm thickness, 80 mm stroke height adjustable backrest, padded in 40 kg/m³ density fire retardant polyurethane foam. With or without 70 mm stroke seat slide.

Armrests:



BR020.N - black nylon height adjustable armrests with PU armpad, 90 mm stroke.

BR021.N - black nylon 3D adjustable armrests with PU armpad, 90 mm stroke.

Mechanism:



MS11H – sincronizzato di colore nero con 5 posizioni di bloccaggio, antishock e pomolo di regolazione intensità oscillazione schienale.

Base:



BG018.N – Ø 700 mm black nylon base with metal ring.



BG004.L – Ø 700 mm polished aluminum base.

Castors:




R26 – Ø 65 mm self-breaking castors for soft floors.



R27 – Ø 65 mm self-breaking castors for hard floors.


Seat slide adjustment.

Pull the lever on the left hand side of the mechanism to adjust the seat depth.



Seat height adjustment.

Pull the right hand side lever while seated to lower the seat height. Pull the lever while standing to increase the seat height.



Tilt adjustment.

Twist the left hand side knob to release the synchronized action. Twist the knob back to lock the synchronized in the desired position.



Tilt tension adjustment.

Twist the right hand side knob to adjust the tilt tension.



Lumbar support adjustment.

Lift or lower the lumbar support in the desired position.

